

Catalog No.: D501P

Product Name: **TaqProbe Direct-qPCR Kit**

Description: **TaqProbe Direct-qPCR Kit** contains all the reagents needed for quick preparation of genomic DNA and the qPCR master mix for **TaqMan real-time PCR assay**. Any type of tissues can be used for this kit, such as mouse tails and nail tips, and other animal tissues and cells; plant leaf, root or seeds; bacteria, fungi and other samples.

Kit Contents:

Size:	125 rxns
DNA Extraction Solution A	13 ml
DNA Extraction Solution B	1.5 ml
2X TaqProbe qPCR MasterMix-P	1.25ml

Storage: The whole kit can be stored at 4°C for up to three months or at -20°C for long-term.

General Protocol

I. DNA Sample Preparation:

- Place the sample into a PCR tube:
 - For mouse tissues from tail, ear or nail: 1-3 mm in length or diameter.
 - Animal tissues: 1-3mg.
 - Cultured cells: 10µl of cell culture.
 - Plant materials: 1-3mg (approximately the size of a sesame seed)
 - Other samples: similar amount or volume as above.
- Add 100µl of the **DNA Extraction Solution A** into the tube containing the sample.
- Heat the sample for 10 minutes at 95°C. This can be easily done in a PCR machine.
- Take out the sample and add 10µl of the **DNA Extraction Solution B**.
- Mix well with vortex or by shaking the tube a few times
- The sample DNA is now ready for real-time PCR or stored at or below 4°C for future applications.

Note:

- The sample can be centrifuged briefly before use.
- Use only the supernatant for qPCR and avoid any undigested tissue or debris.

II. Real-time PCR:

Prepare the qPCR reaction mixture as follows:

Components	Volume 20µl	Final Concentration
TaqProbe qPCR Master Mix	10.0µl	1x
TaqMan Probe	Variable	100-300nM
Forward Primer	Variable	100-500nM
Reverse Primer	Variable	100-500nM
Sample DNA	1.0µl	<500ng
RNase-free Water	Up to 20µl	-
Total Volume	20µl	-

Perform real-time PCR according to your favorable program, or try the following program.

Step	Temperature	Duration – Standard	Duration - Fast	Cycles
Enzyme Activation	95°C	10min	10min	Hold
Denature	95°C	15sec	3sec	40
Anneal/extend	60°C	60sec	30sec	

Recommendations for Optimal Results

- Avoid contamination and avoid repeated freeze-thaw cycles for the reagents.
- TaqProbe qPCR Master Mix components are light sensitive; avoid exposure to light.
- If needed, ROX reference can be used since the qPCR MasterMix-P contains the reference dye ROX.
- Start the real-time PCR as soon as the reaction mixture is prepared and always keep the reaction mixture chilled prior to PCR reactions.

Troubleshooting: Problems and Solutions

Q1. The samples are not completely digested or dissolved?

A1. Samples are not expected to be digested or dissolved completely. Do not worry. Sufficient DNA will be released for PCR without complete digestion of the samples.

Q2. There is little or no real-time PCR signal detected?

A2. Please consider one of the following:

- a) Make sure that there are no PCR components missed.
- b) More PCR cycles may be needed.
- c) Primers may not be designed optimally.
- d) Adjust the real-time PCR parameters to find out the optimal condition for your primers.
- e) Too much sample may have been used. In that case, the samples can be easily diluted 10 times with H₂O or 10mM Tris-HCl buffer, pH 8.5.

Q3. A high background/noise signal?

A3. Adjust your annealing temperature or other parameters for your PCR program.

Q4. The negative control shows false positive signal?

A4. Reagents or your samples may be contaminated.

END