



Catalog No.:	G003
Product Name:	Protein Gel SafeStain
Size:	For 1000ml
Description:	<p>Protein Gel SafeStain is a highly sensitive fluorescent stain designed for detecting proteins in polyacrylamide gels. Compared to traditional Coomassie® blue stain,</p> <p>Protein Gel SafeStain has the following advantages:</p> <ul style="list-style-type: none">• Higher sensitivity: Detecting as little as 0.2 ng protein.• Simple and fast: As short as 60 minutes is needed for staining.• Compatibility: Compatible with MS and protein sequencing.• Wide linear detection range. At least three orders of magnitude.
Detecting:	The Protein Gel SafeStain has two excitation wavelength peaks of 300nm and 460nm; and an emission wavelength of 600nm. So the stained proteins can be viewed and pictured using UV transilluminator, blue light transilluminator or a laser scanner.
Storage:	Store under dark at 4°C or room temperature.

Staining Protocol:

Prepare the Protein Gel SafeStain Solution (1X):

1. Use a clean bottle, add and mix 600ml of water, 300ml of methanol, and 100ml of acetic acid.
2. Transfer the entire content of the **Protein Gel SafeStain** from the tube (1ml) to the above solution.
3. Mix well and store this **Protein Gel SafeStain Solution** at room temperature or at 4 °C, protecting from light.

Protein Gel Staining:

Note: The protocol is optimized for mini protein gels with 1mm thickness. For larger or thicker protein gels, the volume of **Protein Gel SafeStain Solution** should be increased accordingly; and longer incubation time will be needed.

1. **Run** protein gel as usual according to your standard protocol.
2. **Stain** the gel with the 1X **Protein Gel SafeStain Solution** (50-80ml may be needed) at room temperature for 60 min with shaking. **Note:** longer staining will not damage the final result, and sometime will be recommended.
3. **Wash** gel with 100ml of **Protein Gel Wash Solution** (prepared before use as a mixture of 60ml of water, 30ml of methanol and 10ml of acetic acid) for 20 min with shaking. **Note:** longer than 20 min wash will not affect the final result, and sometime will be encouraged.
4. **View and taking picture** of the stained gel using a 300nm UV transilluminator, blue light transilluminator or a laser scanner.

This Product is For Research Use Only

For the proper disposal of this product, follow your University and / Or Company's waste disposal Guidelines.