



Catalog No.: C109

Product Name: Tryptone

Size: 500g

Description: High quality peptone, molecular biology certified (MBC). This product is a pancreatic digest of casein. It is widely used for general bacteriological research work with a variety of microorganisms, such as for making various microbiological media and as an additive for other culture media.

Specifications: Typical representative sample of different lots of peptone.

Amino nitrogen: 4.2%
Total nitrogen: 13.13%
Appearance: yellowish powder
Loss on drying: 3.3%
pH(2% solution): 6-7

Storage: Room temperature or below, under dry conditions.

Note: This Product Is For Research Use Only.

<u>Related Products</u>	<u>Catalog No.</u>
• Terrific Broth	C131
• L-Broth	C122
• LB Broth	C120
• Agar	C110
• Yeast Extract	C108