



Catalog No.: C110

Product Name: Agar

Size: 500g

Description: High quality agar, derived from *Gelidium* seaweed. This product is widely used for general bacteriological research for making various microbiological media and other tissue culture media.

Note: This Product Is For Research Use Only.

Specifications: Typical representative sample of different lots of agar.

Appearance:	white to light cream powder
Particle size:	95% passes 60 mesh sieve
Gel strength:	550-900g/cm ² (1.5%)
Gel point:	32-38°C(1.5%)
Melting point:	82-88°C(1.5%)
pH in gel:	6-7(1.5%)
Solubility:	in boiling water

Storage: Room temperature or below, under dry conditions.

For Research Use Only!

Related Products	Catalog No.
• Terrific Broth	C131
• L-Broth	C122
• LB Broth	C120
• 2xYT Agar	C110
• Yeast Extract	C133