

**Catalog No.:** RNA-SC

**Product Name:** RNA Spin Columns

**Size:** 100/set

**DESCRIPTION:** RNA Spin Columns are silica membrane based centrifugal columns. They have been widely used in RNA purifications, such as total RNA and mRNA purification. These spin columns are also widely used for RNA impurity removal or RNA clean-up in applications like RT-PCR and gene expression study experiments. They are the main components of the RNA mini-prep kits. The Spin Columns are compatible with a wide range of buffers from all brand names.

**Specifications:**

<b>Format</b>	RNA binding column and 2.0ml collection tube as one set.
<b>Principle</b>	Selectively RNA binding
<b>Equipment</b>	Microcentrifuge, Vacuum Manifolds
<b>Binding Capacity</b>	20 µg
<b>Matrix / Construction</b>	Silica-based / Polypropylene
<b>DNase/RNase</b>	None detectable

**Storage:** Stored at room temperature